

Unitron accessory portfolio (Effective February 2026)

Accessory	Image	Compatible Devices	Description	Notes	Pairing
Charger RIC S		Moxi S-R	Offering 24 hours of use on a single charge. Charger RIC S features magnetic insertion for easy management.	All chargers come with 3 year warranty. Utilises USB-C charging.	N/A
ChargerGo RIC S		Moxi S-R	Portable charging for those on the go, ChargerGo RIC S holds 3 additional full charges.	All chargers come with 3 year warranty. Utilises USB-C charging.	N/A
Charger ITE		Insera S-R	Offering up to 29 hours of use on a single charge. Charger ITE features a snap and charge mechanism for easy management.	LED light on the charger contacts indicate charging status. All chargers come with 3 year warranty. Utilises USB-C charging.	N/A
Moxi V-RS charger		All Moxi V-RS devices	Offering 24 hours of use on a single charge. The Moxi V-RS charger utilises magnetic contacts for easy insertion and removal. It's sleek and compact enough to fit into your pocket. Who is it best for? Those who prioritise design and style without having to forego the convenience of rechargeable technology.	All chargers come with a 3 year warranty. Utilises USB-C charging.	N/A
Moxi V charger		Moxi V-R and Moxi V-RT devices	Offering 24 hours of use on a single charge for the confidence to go all day ¹ . The Moxi V charger is equipped with magnetic contacts for easy insertion and removal. Who is it best for? Those who need the easiest and most convenient solution to overcome management and dexterity issues.	All chargers come with a 3 year warranty. Utilises USB-C charging.	N/A

¹ Actual battery life depends on several factors including hearing loss, battery age, and processing demands like direct connectivity.

Charger Moxi R 3		All Blu, Discover Next and Discover Moxi rechargeable devices	<p>Offering 24 hours of use on a single charge for the confidence to go all day.¹ The Charger Moxi R 3 can house a drying capsule (available to order from Unitron) to manage moisture in the hearing instruments.</p> <p>Who is it best for?</p> <p>Perfect for those who have difficulty with batteries.</p>	<p>All chargers come with a 3 year warranty.</p> <p>Utilises Micro-USB charging.</p>	N/A
Charger Combi BTE 2		Stride V-PR, Stride B-PR and DX Stride P R devices	<p>Offering 24 hours of use on a single charge, the Charger Combi BTE 2 features an in-built voltage booster to ensure stable charging.</p> <p>Who is it best for?</p> <p>Perfect for those who have difficulty with batteries.</p>	<p>All chargers come with a 3 year warranty.</p> <p>Utilises Micro-USB charging.</p>	N/A
Power pack		Charger Moxi R 3 and Charger Combi BTE 2	<p>This accessory instantly transforms the Charger Combi BTE 2 or Charger Moxi R 3 into a convenient, portable solution. The Power pack holds 7 full charges, without the need to plug the charging case into an outlet.</p> <p>Who is it best for?</p> <p>For those who require power on the go.</p>	<p>The power pack clips onto the bottom of the Charger Moxi R 3 and Charger Combi BTE 2.</p>	N/A
Stride P R charger		All Vivante, Blu, Discover Next and Discover Stride rechargeable devices	<p>Offering 24 hours of use on a single charge for the confidence to go all day.¹ The Stride P R charger comfortably houses earmoulds.</p> <p>Who is it best for?</p> <p>Perfect for those who have large ear moulds.</p>	<p>All chargers come with a 3 year warranty.</p> <p>Utilises USB-C charging.</p>	N/A
TV Connector		All wireless Smile, Vivante, Blu, Discover Next, Discover and T Moxi All devices.	<p>An easy plug and play solution that works with all TVs, allowing wireless streaming of the TV signal straight to the user's hearing instruments, all the while not interrupting the sound from the TV so users can enjoy watching their favourite shows and movies with others.</p> <p>Who is it best for?</p> <p>Those who want to watch their favourite shows and movies in high-quality stereo sound by wirelessly connecting their hearing instruments to the TV.</p>	<p>Multiple pairs of hearing instruments can connect to the same TV Connector simultaneously.</p> <p>Discover instruments can only be paired to one TV Connector at a time. Blu and Discover Next and Vivante instruments can be paired to two TV Connectors.</p>	<p>Plug the TV Connector into an adequate power source (confirmed by a green light on the LED light indicator).</p> <p>Press the pairing button for 3 seconds and device/s will pair automatically.</p>

PartnerMic		All wireless Smile, Vivante, Blu and Discover Next devices.	<p>The speaker's voice will be transmitted directly into the user's hearing instruments from up to 25 meters away. It uses its 3 microphones to create 2 beams that will focus towards the speech.</p> <p>Who is it best for?</p> <p>The PartnerMic is designed for use in environments where users want to be able to hear a single speaker in loud noise or at a distance.</p>	The PartnerMic must be worn or held by the speaker at all times and positioned within 20cm of the speaker's mouth for optimal performance.	Ensure the PartnerMic is fully charged and switch it on. Press the pairing button for 3 seconds (confirmed by a flashing blue light on the LED light indicator), and device/s will pair automatically.
Remote Control		All wireless Smile, Vivante, Blu and Discover Next devices.	<p>Gives users the ability to easily and discreetly make simple adjustments to volume or programs.</p> <p>Who is it best for?</p> <p>For those who want an easy-to-use remote control or those who cannot or do not want to use the Remote Plus app for adjustments.</p>	The battery is a replaceable cell, not rechargeable. The Remote Control requires one lithium button cell (CR2032).	Hold down the + button and turn the Remote Control on at the same time, this will put the Remote Control into pairing mode (confirmed by a flashing blue light on the LED light indicator), device/s will pair automatically.
IIC remote		Non-wireless Smile, Blu and Discover Next Insera 10A Omni devices.	<p>The pocket-sized IIC remote gives you a choice between magnetically adjusting between programs, volume up only or volume down only (programmable in Truefit).</p> <p>Who is it best for?</p> <p>For those who like the added convenience of having further control over their hearing.</p>	IIC Remote function must be requested at the time of hearing instrument order.	N/A
Remote Plus App	 Download on the App Store GET IT ON Google Play	All wireless Smile, Vivante, Blu, Discover Next and Discover devices.	<p>Users can access the Remote Plus app to make adjustments to essential features and share their hearing experiences with their clinician in real time, it also enhances their hearing experience with on-demand tips, troubleshooting and videos on how to care for and use their hearing instruments.</p> <p>Who is it best for?</p> <p>Those who want to control their hearing instruments and share in-the-moment impressions of their hearing instrument's performance, all from their smartphone.</p>	Requires a minimum of iOS 15 or Android 8 or newer, with Bluetooth 5.0 or newer. Remote Plus app functionality will differ depending on what device is paired.	N/A