


Put the power to
choose in the palm
of their hands

The Remote Plus app delivers a truly personalised experience to every client by empowering them with convenient customisation options.

Remote Plus app



Simplicity + Personalisation

The Remote Plus app provides a modern, user-centric experience that features:


- A vibrant, user-friendly design that encourages client engagement
- Intuitive icons and compelling pictures for the app's features and programs
- Seamless and efficient navigation for in-the-moment personalisation
- Smooth screen flows and enhanced accessibility features

Your clients can choose how they want to hear the world, plus how they like to interact with the app using the following personalisation options:

- Comfort and Clarity optimisation, volume adjustments and program selection all from the Home screen
- Simple adjustments of situation-specific programs and features for unique circumstances and preferences
- The ability to customise program names and choose between program photos and icons
- Two app modes to choose from: Advanced and Classic




Advanced mode



Classic mode


Using the innovative tools found in the Remote Plus app, you can provide more personalised hearing care to your clients along every step of their journey from fitting through ongoing support.



Lifestyle data


Clients get a clear picture of their listening lifestyle allowing them to make decisions with confidence. The Remote Plus app also features an activity tracker and step counter.*

*Availability may be limited by product style.




Connected care

Ratings allows clients to provide in-the-moment feedback, while capturing hearing instrument settings to help make your counselling more efficient and personalised. With **Success Check**, you have access to insightful listening lifestyle data, and you can send fine-tuning adjustments directly to the client using **Remote Adjust**.



Find my Devices

Give your clients peace of mind knowing there's help to locate lost or misplaced hearing instruments. Find my Devices shows the known location of their hearing instruments, and provides guidance based on signal strength, to pinpoint their location.



Everyday support

Send your clients home with confidence knowing they have everything they need right at their fingertips. The app includes support resources including how-to videos, tips, and frequently-asked questions. Your clients receive guidance for the day-to-day use and maintenance of their hearing instruments through **Coach**, the app's virtual hearing aid guide. Messages are delivered based on their individual journey and are accessible anytime by tapping on Coach in the Support tab.



Download the app and try it for yourself using demo mode:



*For the use of the Remote Plus app, Unitron hearing aids with Bluetooth® connectivity are required. The Remote Plus app works on Apple and Android smartphones. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Google Play and the Google Play logo are trademarks of Google Inc. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of each marks by Sonova AG is under license.



Unitron makes life vibrant with amazing hearing solutions designed to make the experience easy. Because everyone deserves to **Love the experience.™**

unitron.com/au
1800 212 313

© 2025 Sonova AG, or its affiliates.
028-6794-48/V3.00/2025-02/jb

