

A new experience awaits

Five steps to a totally personalised hearing solution

Product Selection Guide 4.0

unitron™ Love the
experience

A Sonova brand

Five steps to the hearing solution that's right for you

Step one

Where do you experience hearing challenges?

Step two

Explore the technology levels

Step three

Pick a style that suits you

Step four


Colour your choice

Step five

Connect to convenience

Step one

Where do you experience hearing challenges?






















Step two

Explore the technology levels

Now we'll select the technology level that best matches your listening lifestyle.

	What's included	What you're missing
Level 9 ★★★★★	Premium performance in 7 environments and the most sophisticated features that let you focus on speech and know where it's coming from in crowds and background noise (where you spend ____% of your time). Plus enjoy improved awareness of where sounds come from and the ability to hear in quiet as well as delightful music, automatically.	Nothing. You've got the whole package!
Level 7 ★★★★	Advanced performance in 6 environments with a great feature that will help you focus on speech coming from the front and sides in challenging situations, improved awareness of where sounds come from and the ability to hear in quiet situations, plus an automatic music program.	The most sophisticated feature to focus in different directions in noisy conversations plus automatic adjustment for comfort in very high noise.
Level 5 ★★★	Standard performance in 4 environments that includes a feature to help you to know where sounds are coming from and the ability to hear in quiet environments.	An advanced noise reduction feature to further enhance speech from the front and automatic engagement of a music program.
Level 3 ★★	The essentials in 2 environments. Limited ability to adjust automatically to accommodate multiple, complex listening environments.	Literally everything else.

Improve quiet and challenging conversations with our latest advanced technology.


Level 9 ★ ★ ★ ★ ★	Level 7 ★ ★ ★ ★ ★	Level 5 ★ ★ ★	Level 3 ★ ★
7 environments	6 environments	4 environments	2 environments
 Noise	 Conversation in noise	 Quiet	 Conversation in a crowd
 Conversation in noise	 Music	 Conversation in a small group	 Conversation in quiet
 Music	 Quiet	 Conversation in a crowd	
 Quiet	 Conversation in a small group	 Conversation in quiet	
 Conversation in a small group	 Conversation in a crowd		
 Conversation in a crowd	 Conversation in quiet		
 Conversation in quiet			

Sound performance by technology level

Level 9

Premium

★★★★★




Spatial Awareness

Binaural microphone mode designed to accurately preserve important spatial cues typically lost with hearing aids.

Level 7

Advanced

★★★★




Spatial Awareness

Binaural microphone mode designed to accurately preserve important spatial cues typically lost with hearing aids.

Level 5

Standard

★★★



Spatial Awareness

Binaural microphone mode designed to accurately preserve important spatial cues typically lost with hearing aids.

Level 3


Essential

★★

Pinna Effect


Minimises the impact of acoustic cues typically lost with hearing aids.

High complexity




Speech from front

AutoFocus 360 identifies and focuses on speech from the front to provide a greater increase of speech than other technology levels. Dynamic noise reduction is active and adjustable to further reduce surrounding noise.




Speech from side




Speech from back

AutoFocus 360 can identify and focus on speech from any direction while reducing interfering noise.



Speech from front

AutoFocus identifies, focuses on and increases speech from the front. Dynamic noise reduction is active and adjustable to further reduce surrounding noise.



Speech from side and back

AutoFocus applies a wide front beam to capture speech from the front and partially from the sides.



All high complexity situations

Applies a wide front beam capturing the speech from the front.



All high complexity situations

Applies a wide front beam capturing speech from the front.

Step three

Pick a style that suits you

Choose from products available
on the **Blu** platform.

All of our hearing aids on the Blu platform are stylish, comfortable and easy to use. Most models feature direct connectivity to seamlessly connect via Bluetooth, so you can stream music, podcasts, calls, or digital assistants directly into your hearing aids, controlling your media with just a tap.* And the rechargeable options are quick to charge and proven to last all day.



Moxi B-R

Our smallest rechargeable
with direct connectivity

☐ Receiver-in-canal (RIC)



Moxi B-RT

Our rechargeable with telecoil
and direct connectivity

☐ Receiver-in-canal (RIC)



Moxi B-312

Our smallest RIC with direct connectivity and
traditional 312 disposable battery

☐ Receiver-in-canal (RIC)

Charging case for
Moxi B-R and Moxi B-RT
rechargeable
hearing aids





Stride™ B-PR

Our rechargeable BTE with direct connectivity

☐ Behind-the-ear (BTE)



Stride B-UP

Our Ultra Power BTE with direct connectivity, telecoil and traditional 675 disposable battery

☐ Behind-the-ear (BTE)

Charging case for
Stride rechargeable
hearing aids



Insera™ B-312

Our wireless ITE model with direct connectivity

☐ In-the-ear (ITE)



Insera B-312 NW O

Our non-wireless omni ITE model with a telecoil

☐ In-the-ear (ITE)



Insera B-10 NW O

Our smallest non-wireless omni ITE model with a telecoil

☐ In-the-ear (ITE)

Step four

Colour your choice



Sand
(R8)



Amber
(P2)



Espresso
(P4)



Platinum
(P6)



Pewter
(P7)



Charcoal
(P8)



Beige
(26)



Tan
(14)



Cocoa
(22)



Brown
(28)



Black
(06)*

*Available for Insignia B-312 and Insignia B-10 IIC styles only.
Moxi B-R shown at actual size.

Step five

Connect to convenience



TV Connector

Stream your favourite media wirelessly in stereo sound with the TV Connector.



PartnerMic

PartnerMic is an easy-to-use lapel microphone ideal for one-on-one conversations.



Remote Control

Discreetly adjust volume and switch programs with the Remote Control.

Remote Plus app

The Remote Plus app puts easy personalisation at your fingertips.



Scan to download the Remote Plus app.



iOS



Android

Need more help in challenging listening situations like meetings or conversations in noise or over distance?

Let's discuss if a Roger™ device could be a good fit for you.

Accessories shown at actual size.

Accessories shown are available with wireless products on the Blu platform.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Google Play and the Google Play logo are trademarks of Google Inc

What's next?

- ☐ Discuss your needs
- ☐ Try at home
- ☐ Come in for a follow-up assessment
- ☐ Personalise your solution
- ☐ Come in for your final fitting on:

Recommendation:			
Technology level:	Style:	Warranty:	Cost:
Accessories:			

Notes:



At Unitron, we empower people with life-enhancing hearing experiences that fit seamlessly into their world. Our sound performance technology, Experience Innovations, and intuitive design work perfectly together for unmatched personalisation and optimisation. Because everyone deserves to **Love the experience.**TM

© 2022 Unitron. All rights reserved.

028-6543-02/v4.00/2022-09/dr

unitron.com/au



sonova
HEAR THE WORLD